

INDICATIONS

SPORIGERM FAZZOLETTI⁺

is a ready to use medical device with an effective disinfectant action against Gram-positive and Gram-negative bacteria, yeasts, Mycobacterium Tuberculosis (TBC), Viruses (including HIV-HBV-HCV).

A convenient solution which makes

SPORIGERM FAZZOLETTI⁺

a specific product for the disinfection and cleansing suitable for surfaces of non-invasive medical devices such as dental chairs, worktops, armrests, handles, turbines, face masks, headrests etc., even in the presence of organic material (blood, saliva, exudates).

INSTRUCTIONS FOR USE

Wipe surfaces and objects to be cleaned, decontaminated and disinfected making sure these are completely dampened. Rinsing is not necessary.

CONTACT TIME

Enveloped viruses (HIV-HBV-HCV).....	5 min
Bacteria	5 min
Fungi (yeasts)	5 min
Non-enveloped viruses (Adenovirus)	5 min
Mycobacterium Tuberculosis (TBC)	60 min

PACKAGING

Sporigerm Fazzoletti⁺

Code 501002.....12 flow packs, 50 wipes each

Primary packaging in compliance with the technical specifications established by Farmacopea.


These materials do not contain latex and are fully compatible with all the components of the formulated.

COMPOSITION 100 grams of product

Benzalkonium Chloride.....	0.50 g;
O-Phenyl phenol.....	0.02 g;
Isopropyl alcohol, Ethyl alcohol, Co-formulants;	
Purified water, q.s. at	100.00 g;

Disinfectant detergent
with tuberculocidal, bactericidal,
fungicidal, virucidal action

Medical Device Class IIb
CE n° 0546



INFORMATION FOR THE EXCLUSIVE USE OF HEALTHCARE PROFESSIONALS

Manufactured by:

GERMO S.p.A.

Via Giotto n 19/21 - 20032 Cormano - MI

Phone +39 02 663 01 938

www.germo.it - info@germodis.com



Distributed by:

IDS S.p.A.

Via S. Cristoforo 28/10

17100 Savona - Italy

www.ids dental.it - info@idsdental.it

Disinfectant detergent
with tuberculocidal, bactericidal,
fungicidal, virucidal action

CHARACTERISTICS

Wet wipes

ACTION

The Quaternary Ammonium Salt, together with the O-Phenyl phenol, in alcohol solution, guarantees a full range of action towards Gram-positive and Gram-negative bacteria, yeasts, mycobacteria (TBC) and viruses (Adenovirus, HIV, HBV, HCV).

EN13727 - EN 13697

Bactericidal 5 minuti

EN13624 - EN 13697

Fungicidal 5 minuti

EN14476

Virucidal (Adenovirus, HIV-HCV, HBV) 5 minuti

EN14348

Mycobacterium Tuberculosis (TBC) 60 minuti


PRODUCT VALIDITY

3 years. The validity period indicated refers to a product stored in its package and correctly used and preserved.

STORAGE

To be stored in cool, dry areas, away from heat sources.

PATHOGENIC AGENTS CLASSIFIED ACCORDING TO THEIR RESISTANCE TO DISINFECTANTS



PRECAUTIONS

Not to be used on materials which can be damaged by alcohol.

QUALITY CONTROL

All components (raw materials, containers, labels etc.) together with the processing phases of each batch are punctually and thoroughly checked internally according to the corporate UNI EN ISO 9001 and UNI EN ISO 13485 certified Quality Management System

AUTHORISATIONS

Class IIa CE 0546 medical device complying with Directive 93/42/EEC as amended by Directive 2007/47/EE

INFORMATION FOR THE EXCLUSIVE USE OF HEALTHCARE PROFESSIONALS

Manufactured by:

GERMO S.p.A.

Via Giotto n 19/21 - 20032 Cormano - MI

Phone +39 02 663 01 938

www.germo.it - info@germodis.com

Distributed by:

IDS S.p.A.

Via S. Cristoforo 28/10

17100 Savona - Italy

www.ids dental.it - info@idsdental.it