

INDICATIONS

SPORIGERM PERACETIC PLUS

is a disinfecting/sterilising powder for invasive surgical instruments.

The active ingredient, the peracetic acid, is formed "in situ" after having dissolved the powder in water to obtain an aqueous solution of peracetic acid. SPORIGERM PERACETIC PLUS can also be used on endoscope reprocessors and is effective for the fast disinfection and sterilisation of: endoscopes, optic fibre, surgical and dental instruments, transducers.

INSTRUCTIONS FOR USE

Dilute SPORIGERM PERACETIC PLUS using lukewarm fresh water at a temperature of approx. 35° C; the solution thus obtained will appear clear, light blue in colour and ready for use after 15 minutes. The solution, once completely colourless, can be used to immerse tools previously washed in the following way:

To decontaminate/disinfect:

Use 8 g of powder (equal to one dosing spoon) for each litre of water for 10 minutes. After contact time, remove and rinse thoroughly.

To sterilise:

Use 16 g of powder (equal to two dosing spoons) for each litre of water for 10 minutes. Subsequently remove the instruments and rinse in sterile water. For aseptic storage use suitable product. The solution is active for a maximum time of 24 hours from preparation. After use, the diluted product can be disposed of as non-hazardous waste.

PACKAGING

Sporigerm Peracetic Plus

Code 502002.....2 jars, 500 g each


Primary packaging in compliance with the technical specifications established by Farmacopea. These materials do not contain latex and are fully compatible with all the components of the formulated.

COMPOSITION

100 g of product
Tetraacetylenediamine..... 25.00 g;
Sodium percarbonate 42.00 g;
Co-formulants q.s. at 100.00 g;

Cold steriliser disinfectant
**active in 10 minutes on spores,
HIV, HCV, HBV, TBC, fungi and bacteria**

Medical Device Class IIb
CE n° 0546



INFORMATION FOR THE EXCLUSIVE USE OF HEALTHCARE PROFESSIONALS

Manufactured by:
GERMO S.p.A.
Via Giotto n 19/21 - 20032 Cormano - MI
Tel. 02 663 01 938
www.germo.it - info@germodis.com



Distributed by:
IDS S.p.A.
Via S. Cristoforo 28/10
17100 Savona - Italy
www.ids dental.it - info@idsdental.it

Cold steriliser disinfectant

**active in 10 minutes on spores,
HIV, HCV, HBV, TBC, fungi and bacteria**

PATHOGENIC AGENTS CLASSIFIED ACCORDING TO THEIR RESISTANCE TO DISINFECTANTS

CHARACTERISTICS

Granular powder white in colour, might have light blue shades.

ACTION

Results of the antimicrobial effectiveness tests are summarised in the following table:

EN13727; EN14561

Bactericidal0.4%..... 10 min

EN13624; EN14562

Yeasticidal (C. Albicans)0.4%..... 10 min

EN14563

Mycobactericidal/ Tuberculocidal 0.4% 10 min

EN13704

Sporicidal0.8%..... 10 min

EN13624; EN14562

Fungicidal (A. brasiliensis) 1.6%..... 10 min

EN14476


Virucidal1.6%..... 10 min

EN14937

Cold steriliser1.6%..... 10 min

PRODUCT VALIDITY

The product has a validity of 2 years from the manufacturing date, it is to be stored in its original package. The product is affected by moisture



STORAGE

To be stored in cool, dry areas away from heat sources.

QUALITY CONTROL

All components (raw materials, containers, labels etc.) together with the processing phases of each batch are punctually and thoroughly checked internally according to the corporate UNI EN ISO 9001 and UNI EN ISO 13485 certified Quality Management System.

AUTHORISATIONS

Class IIb CE 0546 medical device complying with Directive 93/42/EEC as amended by directive 2007/47/EEC

INFORMATION FOR THE EXCLUSIVE USE OF HEALTHCARE PROFESSIONALS

Manufactured by:

GERMO S.p.A.

Via Giotto n 19/21 - 20032 Cormano - MI

Tel. 02 663 01 938

www.germo.it - info@germodis.com

Distributed by:

IDS S.p.A.

Via S. Cristoforo 28/10

17100 Savona - Italy

www.ids dental.it - info@idsdental.it