

Highly effective surface disinfectant spray designed for use on non-invasive medical device surfaces, including those that are sensitive to alcohol. Contains detergent for added cleaning power. Effective against Bacteria, Yeast, Fungi and Viruses including: HBV, HIV and Rotavirus*

Medical Device Class IIb
CE n° 1639

- The original cartridge based delivery system
- Guaranteed and accurate dosing every time
- Safest system available
- Contact time from only 30 seconds
- Wide effective spectrum
- Wide material compatibility
- Contains detergent

KEY FEATURES

Viroxid® Spray is a concentrated, cartridge based delivery system for alcohol free surface disinfectant.

Designed to remove the need to measure and dispense dangerous chemicals and ultimately save practices both time and money while ensuring that the correct exact dosage is used each and every time - improving safety for both staff and patients.

Each Viroxid® Spray cartridge contains exactly the right amount of concentrate for one 500ml spray bottle, just add water.


The cartridge is placed inside the neck of the exclusively designed bottle and is automatically activated by screwing the trigger mechanism into place.



Viroxid® Spray can be used on almost all surfaces within the dental surgery, including non invasive medical devices. Each cartridge dispenser pack contains 10 x 8ml cartridges which is the equivalent of a 5L/5kg container of ready to use disinfectant. Packs can easily be wall mounted with the included sticky adhesive pads. The outer surface of the cartridge dispenser pack is coated in Biocote - a unique antimicrobial, silver ion coating, ensuring the dispenser pack can be wall mounted even in designated 'clean' areas.

COMPOSITION

Aqua, Benzyl-C12-18-Alkyldimethyl Chlorides, C12-14-Alkyl[(Ethylphenyl)Methyl]Dimethyl Chlorides, Didecylmethylammonium Chloride, 2-Propanol (not an active, <1%), Alcohols (not an active) C9-11 Ethoxylated, Dipropylene glycol monomethyl ether, Polyhexamethylene biguanide hydrochloride, Cocamidopropyl Betaine, Tetrasodium EDTA, Dye.



INFORMATION FOR THE EXCLUSIVE USE OF HEALTHCARE PROFESSIONALS

Manufactured by:
Topdental (Products) Ltd.
12 Ryefield Way, Silsden, West
Yorkshire. BD20 0EF. England UK.
Tel +44 (0)1535 652 750 www.topdental.co.uk



Distributed exclusively by:
IDS S.p.A.
Via S. Cristoforo 28/10
17100 Savona - Italy
www.ids dental.it - info@ids dental.it

**EFFECTIVE SPECTRUM
AND CONTACT TIMES**

MICROORGANISMO	TEMPO
Bactericidal	
Pseudomonas aeruginosa	30 sec
Stafilococco aureo	30 sec
Enterococco hirae	30 sec
VRE	1 min
MRSA	1 min
Legionella pneumophila	5 mins
Yeasticidal	
Candida Albicans	30 sec
Enveloped Virus	
Vaccinia virus	2 min
Poxviridae	2 min
Herpesviridae	2 min
Filoviridae (e.g. Ebola, Marburg)	2 min
Flavivirus	2 min
Hepatitis B virus (HBV)	2 min
Hepatitis C virus (HCV)	2 min
Hepatitis Delta virus (HDV)	2 min
Influenza virus	2 min
Paramyxoviridae	2 min
Rubella virus	2 min
Measles virus	2 min
Rabies virus	2 min
Coronavirus	2 min
Human immunodeficiency virus (HIV)	2 min
Human T cell leukemia virus (HTLV)	2 min
Fungicidal	
Aspergillus niger	30 sec
Mycobacteria	
Mycobacterium avium	2 mins
Spores	
Clostridium difficile (C-diff)	5 mins


*Full effective spectrum available at www.topdental.co.uk

TEST STANDARDS

The following test standards have been achieved:

- EN 14885,
- EN 14476,
- EN 14348,
- EN 13727,
- EN 13704,
- EN 13624.

**PATHOGENIC AGENTS CLASSIFIED ACCORDING TO
THEIR RESISTANCE TO DISINFECTANTS**



INSTRUCTIONS FOR USE

Fill the system bottle with a maximum of 500ml clean cold water.
Put 1 a cartridge into the neck of the system bottle.
Place trigger spray onto system neck by passing the dip tube through the Viroxid®
Spray cartridge and screw down.
Gently shake the bottle for a few seconds prior to use.
Viroxid® Spray will remain active for 12 months following activation.
Spray surfaces and objects to be disinfected from a distance of approximately 30 cm so that they are completely moistened and let the disinfectant take action over a period of 30 seconds (to achieve full efficacy leave for 5 minutes).
Allow surface spray to dry or wipe with clean tissue. For professional use only.

PACKAGING

Viroxid® Spray
Code 070100 (10) cartridges dispenser 8 ml (5 liters)
(1) empty 500ml bottle with trigger

Code 070101 (10) cartridges dispenser 8 ml (5 liters)

All primary packaging are compliant to the technical specifications provided by the Pharmacopoeia.
The material used does not contain latex and it is perfectly compatible with all the components of the formulation.

**INFORMATION FOR THE EXCLUSIVE USE
OF HEALTHCARE PROFESSIONALS**

Manufactured by:
Topdental (Products) Ltd.
12 Ryefield Way, Silsden, West
Yorkshire. BD20 0EF. England UK.
Tel +44 (0)1535 652 750 www.topdental.co.uk

Distributed exclusively by:
IDS S.p.A.
Via S. Cristoforo 28/10
17100 Savona - Italy
www.ids dental.it - info@ids dental.it