# I AQUA PRINT MONO (fast set)

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# GB AQUA PRINT MONO (fast set)

Hydrophilic Fast Setting One Phase Impression Material

Instructions for Use

AQUA PRINT MONO (fast set) is a high hydrophilic and fast setting one phase impression material for crowns, bridges, inlays, partial and complete dentures and denture repairs. It can also be used for impression materials for relining, functional impressions, impressions for prosthetic restorations and for impressions using the copper band technique.

AQUA PRINT MONO is a medium viscosity, thixotropic and hydrophilic precision impression material based on addition curing vinyl silicones. The mixing ratio is 5 volumes base paste and 1 volume catalyst paste. The material shows high hydrophilic properties and is compatible for moist surfaces. Therefore it yields an exact reproduction of subgingival areas and other moist surfaces. AQUA PRINT MONO is delivered in

MAXIMIX-cartridges.

The cartridges can be used with suitable mixing machines.

AQUA PRINT MONO is thixotropic. That means it stays in position in the tray and under pressure it flows. The material meets the requirements of ISO 4823 (2000), type 2 (medium-bodied consistency).

## Preparation of the cartridge

Use only original MIXPAC™ mixing cannulas type MBD 381-05-HS. Hold the cartridge with the sealed outlet openings facing upwards. Take the handle of the sealing cap and remove the sealing cap upwards. Avoid excessive force. Insert the cartridge into the mixing machine (according to manufacturer instructions). Attach a new mixing tip to the cartridge. After the mixing tip

is correctly seated (the central hexagonal socket at the lower end of the tip must be properly aligned with the hexagonal drift shaft), place the yellow Bayonet ring over the mixing tip and push it down completely. Turn the ring clockwise (1/4 turn) to lock the mixing tip securely in place. Check whether the hexagonal drift shaft does engage when starting. The spiral in the mixing tip must rotate when starting.

**If the used mixing machine offers different mixing speed levels, use always the lowest mixing speed.**

### Notes:

Discard the first 3 cm of impression material extruded from newly filled cartridges.

Store used cartridge with fixed used mixing tip. It serves as a cap.

To remove the used mixing tip turn the Bayonet ring anti-clockwise and draw it off.

### Impression taking

Load the AQUA PRINT MONO onto the tray. Insert the tray not later than 1:30 minutes (working time at 23°C or 74°F) after extruding the material into the mouth. Maintain position for at least 2 minutes.

### Casting

The impression should be stored at normal room temperature (max. 25 °C). The impression can be poured 30 minutes after removal from the mouth. No other time limitation does exist.

### Electroplating

All impressions can be copper- or silverplated.

### Desinfection

Impressions can be desinfected with standard desinfection solutions.

### Storage

Do not store above 25 °C (78 °F). Do not use after expiration date.

#### Additional notes:

- Do not use with standard (condensation reaction type) silicone impression materials.
- Do not use Latex gloves, they may interfere with the setting reaction of the impression material.
- MIXPAC™ is a registered trade mark of Sulzer Mixpac AG Switzerland.

#### Warranty

IDS Dental warrants this product will be free from defects in material and

manufacture. IDS Dental makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and IDS Dental's sole obligation shall be repair or replacement of the IDS Dental.

#### Limitation of Liability

Except where prohibited by law, IDS will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

**Keep away from children!  
For dental use only!**

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- MIXPAC™ ist ein eingetragenes Warenzeichen der Sulzer Mixpac AG Switzerland.