I AQUA PRINT REGULAR BODY (AUTOMIX, fast set)

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GB AQUA PRINT REGULAR BODY (AUTOMIX, fast set)

Hydrophilic fast setting material for the (one step)

two phase impression technique as well as one phase impression material

Instructions for Use
Rev1 - 17 Gen 2014

AQUA PRINT REGULAR BODY is very versatile. It can be used:

- as a **hydrophilic syringable material** for the (one step) two phase impression technique in conjunction with AQUA PRINT PUTTY
- as a **Correction material** in the (two step) putty/wash technique
- as a syringable **One Phase Impression Material** for crowns, bridges, inlays, partial and complete dentures and denture repairs. It also serves as an impression material for relining, functional impressions, impressions for prothetic restorations and for

impressions using the copper band technique.

AQUA PRINT REGULAR BODY is a **light to medium viscosity hydrophilic** precision impression material based on addition curing vinyl silicones delivered in **AUTOMIX-cartridges**. It can be extruded directly onto the teeth or the tray. AQUA PRINT REGULAR BODY shows **strong hydrophilic** properties. Thus improves the wettability and yields an exact reproduction of subgingival areas and other moist surfaces.

AQUA PRINT REGULAR BODY is thixotropic. That means it stays in the position in the tray, but under

pressure it flows. It meets the requirements of **ISO 4823** (2000) **Type 3** (light-bodied consistency).

Gun load and extruding

Load the gun with the AUTOMIX CARTRIDGE and mixing cannula. For direct extruding onto teeth use "intra oral tips".

Impression taking

Load the auto mixed material onto the tray or directly onto the teeth. Insert the tray not later than **1:30 min** (working time at 23°C or 74°F) after extruding AQUA PRINT REGULAR BODY into the mouth. Maintain position for at least **2 minutes**.

Notes:

Store used cartridge with fixed used mixing cannula. **Do not re-use the cartridge cap.** Closing the cartridge by turning on the cap between use could cause cross contamination of components and thus premature setting.

Casting

The impression should be stored under normal room conditions (max. 25 °C). The impression can be poured **30 minutes** after removal from the mouth. No other time limitation does exist.

Electroplating

All impressions can be copper- or silverplated.

Desinfection

Impressions can be desinfected with standard desinfection solutions.

Storage

Do not store above 25 °C (78 °F). Do not use after expiry date.

Additional notes:

Do not use with standard (condensation reaction type) silicone impression materials.

Do not use Latex gloves, they may interfere with the setting reaction of the impression material.

Warranty

IDS Dental warrants this product will be free from defects in material and manufacture. IDS makes no other warranties including any implied warranty of merchantability

or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and IDS's sole obligation shall be repair or replacement of the IDS product.

Limitation of Liability

Except where prohibited by law, IDS will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

**Keep away from children!
For dental use only!**

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