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## GB Pre operatory gel for dental hygiene and periodontal pocket cleansing

INSTRUCTIONS FOR USE

Read all this leaflet carefully before you start using this product. Keep this leaflet, you may need to read it again.

### 1. WHAT IS CLEAN&WHITE OXYGEL AND WHEN USE IT

Clean&White OxyGel is a pre operatory gel based on 11% Carbamide Peroxide and 1% Disodium EDTA to be used sub and supra gingivally to dissolve plaque and soften tartar deposits.

Before starting the scaling procedures, the gel is easily flowed inside the sulcus by just placing the tip on the gum border and injecting the gel. In few seconds the gel will start releasing Hydrogen Peroxide which drastically controls the bio-load of the sulcus, while the combined action of Urea and Disodium EDTA will produce soluble complexes with calcium and heavy metal facilitating the dissolution of microbial plaque and tartar. This pretreatment will allow speed-up a lot the hygiene procedure working both with hand instrumentation and with no or minimal water spray during ultrasonic scaling, thus dramatically reducing the production of aerosols.

Hydrogen peroxide has a mild hemostatic action. Clean&White OxyGel formula also includes Sodium Hyaluronate to improve the healing of gums after prophylaxis and Potassium Nitrate to reduce the sensitivity from the opening of dental tubules during the scaling. Dental hygiene session will be much faster and comfortable for the patient.

Clean&White OxyGel can be used to brush teeth surface with a prophy brush or cup to remove effectively biofilm and prepare the tooth surface. This will speed up stain removal with air prophy and improve active oxygen penetration during in office teeth whitening procedures.

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