BLANCCONE

Clean&White Xpower

Ultrasonic Scaler & Air Polisher

Instruction Manual **EN** **3**

Manuale d'istruzioni **IT** **9**

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Manuel d'instructions **FR** **21**

Manual de instrucciones **ES** **27**

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Other languages

Video Tutorial

BLANCONE[®]
Clean & White
Xpower

Ultrasonic Scaler & Air Polisher

BlancOne® C&W XPOWER

CE 0123

Precautions

Warning: If you neglect these precautions, it might cause personal injury such as electric shock, fire or damage to the product.

- This product is not a home appliance equipment. It is only applicable for hospitals and dental clinics. The use of the product must comply with the requirements of the relevant operating specifications and relevant laws and regulations of the medical department. The user must be a professionally trained and qualified dentist or technician. Adequate protection (such as goggles, mask, etc.) should be equipped during operation to prevent cross-infection.
- Please use an independent power outlet with protective grounding. Do not unplug the power cord with wet hands.
- Please fully insert the power plug of the device into the power socket, please do not use other power source other than the specified voltage.
- Do not damage, modify, pull, bend or twist the power cord excessively, and do not place heavy objects on the power cord.
- After using the machine, turn off the power switch and pull out the power plug.
- When the equipment recovers after power off, please await the equipment to be stable before stepping on the pedal. The operator should release the pedal in time when the power is off.
- Please do not place the product on an unstable workbench, such as a swaying table, an inclined surface, or a location subject to vibration.
- This equipment is a reusable product. Keep the equipment clean before and after use. Before each use, the handpieces, tips, torque wrench etc. must be disinfected and sterilized. It is suggested to disinfect and sterilize according to the recommended method in Chapter 5 of this manual.
- Do not knock or scrape the handpiece. Do not pull the cable forcefully during the use of the equipment to avoid damages.
- Do not bend or polish the tips. The tips must be screwed to the handpiece with a torque wrench, and water spray must be generated during operation.
- Before each operation, you should work outside the patient's mouth for more than 10 seconds under the condition of water to drain the water in the pipeline at the back of the handpiece
- When the tip is damaged or worn out, the vibration intensity might decrease. The operator should replace a new tip in time according to the clinical situation. It is not recommended to use the tip after the abrasion of tip exceed 2mm or more. You can use the tip comparison indicator to check. If the length is out of the first line, the power remains the same; if it is between the first and second line, the power would decreased to 80%; if it reaches or surpass the second line, the power would drop massively, on the occasion, replacement of tip is recommended.
- Do not screw or unscrew the tip while stepping on the foot switch, or when the machine is on working.
- The inner thread of the tip manufactured by certain manufacturers are rough, rusty and might break the thread system, which will cause irreparable damage to the scaler. Please use corresponding tips of Clean & White brand.
- Choose the appropriate power according to different types of tips. See Annex Tip power chart There are recommended power and water volume for different types of tip in the system. You can double-click the current mode to choose the tip.
- Do not use unclean water
- Improper cleaning and treatment of titanium implants, porcelain restorations, etc. can easily cause loosening of the adhesive, cracking of the porcelain restorations, or even cracking of the porcelain. The cleaning or treatment of the oral cavity should be carefully considered.
- Noise will be generated during the scaling process. Those who are sensitive to noise can wear earplugs by themselves.
- We will be responsible for the safety only when the maintenance, repair and modification of the machine are carried out by the company or its authorized distributors, the replacement parts are from the company, and the operation is on the basis of instruction manual.
- This equipment has electromagnetic interference. Do not use it around electronic surgery, and be cautious in an environment with strong electromagnetic interference when using this equipment.
- This equipment does not contain toxic or hazardous substances, and it can be discarded in accordance with the relevant laws and regulations on discarded medical devices.
- Do not step on the foot pedal when the cable of the air polish handpiece is removed from the main unit.
- Before using the ultrasonic handpiece, please make sure that the cable of the air polish handpiece is correctly placed on the handle bracket; similarly, before using the air polish handpiece, please make sure that the cable of the ultrasonic handpiece is correctly placed on the ultrasonic handle bracket.
- The tip must be tightened.
- The air polishing handpiece nozzle should not be aimed at people.
- There might be injury if the powder accidentally sprays into the eyes. We strongly recommend that all personnel (doctors, nurses, patients) wear goggles during air polishing treatment.
- Please do not unload the powder tank when stepping on the pedal or when the machine is working.
- Before replacing the air polishing handpiece or the nozzle, please use a syringe to blow the moisture at the joints at both ends (especially the gas interface) to prevent moisture from entering the gas path and avoid clogging of the powder in the pipeline.
- The device must not be used in MRI environment for the device is easily affected by the electromagnetic emission and would not work or work normally.
- When you meet circumstances where the main unit expires, the misuses lead to the short circuit of the circuit board or accidentally dropping the device results in the damage of components, the device should no longer be reused.

Contraindications

- The hemophilia patient is forbidden to use this equipment
- The patients or the doctors with heart pacemaker are forbidden to use this equipment.

- Heart disease patients, pregnant women and children should be cautious to use the equipment.
- Patients with respiratory diseases such as asthma and chronic bronchitis are not allowed to use this device.

Symbol Instruction

Symbol	Implication	Symbol	Implication
	Trademark		Temperature limit: -20°C - +55°C
	Consult instructions for use		Atmospheric pressure limitation: 70kPa-106kPa
	Protective earth (ground)		Humidity limitation: 10%-93%
	Caution		Keep dry
IPX1	Anti-drip equipment		Fragile, handle with care
	Sterilizable in a steam sterilizer (autoclave) at 134°C		Waste electrical and electronic equipment
	Recovery	BOOST	Boost mode
	Manufacturer		Anhydrous mode
	For indoor use only		Clean mode
	Date of manufacture	CE 0123	CE marking with identification number of the Notified Body
	Type B applied part	EC REP	Authorized representative in the European Community
	Alternating current	SN	Serial number
	ON (power)	MD	Medical device Authorized representative in the United Kingdom
	OFF (power)	LOT	Batch number
ON/OFF	ON/OFF	UDI	Unique device identifier
	Functional earth	UK RP	Authorized representative in the United Kingdom
	Foot switch		

Forward

Guilin Refine Medical Instrument Co., Ltd. is a high-tech enterprise integrating R&D, production and sales of dental products, with a complete quality assurance system.

1 Product introduction

1.1 Brief introduction

Xpower Ultrasonic Scaler & Air Polisher combines ultrasound scaler and air polishing scaler. The main features are:

- According to the selected handpiece, automatically switch the working mode.
- Touch LCD screen.
- Special chemical solutions such as hydrogen peroxide, sodium hypochlorite and chlorhexidine can be used to improve the clinical treatment effect.
- LED handpiece enables convenient Clinical operation.
- The automatic frequency tracking system for the best working condition, which brings more stable performance of device.
- Detachable ultrasonic handpiece and air polishing handpiece can be sterilized under high temperature of 134 and high pressure of 0.22MPa.

1.2 Model: Xpower

1.3 Intended use

This product includes an ultrasonic system and an air polish system. The ultrasonic system is used for periodontal treatment in oral clinical treatment to remove tartar and plaque, root canals cleaning; Air polish system is used to remove plaque and pigment, it can also used for implants maintenance.

1.4 Device configuration

Please refer to the packing list for device configurations.

1.5 Structure and components

The machine consists of main unit, handpiece (ultrasonic handpiece, air polishing handpiece), cable, bottle, dental air polishing powder, torque wrench and foot pedal, Power cord.

1.6 Replacement instruction for main accessories.

Name	Model	Replacement instruction
Air polishing handpiece	PA / PAP	Reusable and Durable; Please change a new handpiece when the handpiece is blocked and cannot be dredged, or the sealing performance is reduced, or the ventilation pipe is seriously worn out
Ultrasonic handpiece	S	Reusable and Durable; Please change a new handpiece when the output power drops significantly or does not work.
Power cord	/	Reusable and Durable; Please replace the power cord when it is cracked or non-conductive
Water pipe	/	Reusable and Durable; Please change a new pipe when the water pipe is cracked and leaked.
Torque wrench	/	Reusable and Durable. Please change a new one as per Instruction Manual when its broken.
Foot pedal	RFS05/RFS02	Reusable and Durable. Please change a new one as per Instruction Manual when its broken.
Water bottle	/	Reusable and Durable. Please change a new one as per Instruction Manual when its broken.
Powder tank	/	Reusable and Durable. Please change a new one as per Instruction Manual when its broken.

1.7 Device safety classification

- 1) Classified by operation mode: Continuous operating device
- 2) Type of protection against electric shock: Class I
- 3) Degree of protection against electric shock: B type
- 4) Applied part of the equipment: tip
- 5) Degree of protection against harmful ingress of water: Ordinary equipment (IPX0). Foot pedal is anti-drip device (IPX1)
- 6) Degree of safety application in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide:
Equipment cannot be used in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide.

1.8 Main technical specification


Specification	Xpower
Tip thread specifications for compatibility	M3*0.5
Main Unit Size (mm)	320mm*330mm*120mm
Length*Width*Height	
Main unit weight	4.0kg
Handpiece model	S - PA - PAP
Control mode	Touch control
Water supply mode	Water tank
Power adapter input	100VAC - 240VAC 50Hz/60Hz
Input power	200VA
Main unit fuse	T1.6AH 250V
Foot pedal shell protection Level	IIPX1
Tip output features	30kHz±5kHz; Output Tip output features main vibration offset of tip 1µm ~ 90µm; Semi-offset force 0.1-5N; Output power of tip 3W-20W
Water inlet pressure	0.1bar-5bar (0.01MPa-0.5MPa)
Air inlet pressure	5.5bar-7.5bar (0.55MPa-0.75MPa)
Water temperature	≤ 45°C
Software version	1.0.0

1.9 Operation environment


Environmental temperature: +5°C - +40°C ,
Atmospheric pressure: 70kPa - 106kPa,
Relative humidity: 30% - 75%,
Cooling water temperature: +5°C - +25°C

2 Installation

2.1 Main unit diagram



2.2 Touch panel


2.3 Symbol Implication

Symbol	Implication	Symbol	Implication	Symbol	Implication
Tip P1	Tip currently in use	G	Scaling mode	P	Periodontal mode
E	Endo mode	A	Air polishing mode	C	Cleaning mode
Document icon	Data storage function	Gear icon	Setting	Water heating icon	Water heating switch
Water drop icon	Turn down water volume/ power/ air pressure	Water drop icon with plus	Turn up water volume / power / air pressure	Bar chart icon	Power adjustment
Water drop icon	Water volume				

2.4 Schematic diagram of handpiece

2.5 Schematic diagram of tip installation

2.6 Installation procedures

- 2.6.1 Open the package, check whether the equipment is complete as per the packing list, and place the main unit on a solid place, facing the operator.
- 2.6.2 Plug the external air pipe connector into the air intake connector on the back of the main unit (as showed in Figure 8).
- 2.6.3 Fill the water bottle with an appropriate amount of water, and fix the water bottle into the sink seat directly above the main unit (it is recommended to smear a layer of petroleum jelly on the O-ring at the bottom of the water bottle to facilitate the plugging and unplugging of the water bottle)
- 2.6.4 Installation of foot pedal
- 1) Wired mode: Connect one end to the foot pedal port of the foot pedal, and the other end to the foot pedal port on the main unit.
 - 2) Wireless mode: Remove the pedal battery cover and install 2 AA batteries. Switch on the main unit and the foot pedal can be connected and controlled (when leaving the factory, the Bluetooth foot pedal and the main unit have been paired and connected. If you change the foot pedal, please enter the "System Settings" to match the Bluetooth foot pedal.)
- 2.6.5 Connect the ultrasonic handpiece and the air polishing handpiece into the corresponding cables, and place the handpiece on the brackets on both sides of the main unit. The ultrasonic handpiece is on the left and the air polishing one is on the right.
- 2.6.6 Turn off the power switch, connect the power cord to the power socket on the back of the main unit, and then plug in the power, as shown in Figure 2.
- Warning 1:** When the power cord is connected to the network power supply, it must be connected to a protective ground; When the power cord is connected to the network power supply, do not place or install the product where it is difficult to disconnect the network power supply.

3 Function and operation

3.1 Description of the foot pedal function

Different machines are equipped with different kinds of foot pedal, please refer to the actual packing list.

3.1.1 Multifunctional foot pedal

- a) Choose the mode of foot pedal according to 2.6.4, place it in the flat ground
- b) The pedal and its functions are shown below:

Button	Operating mode	Function	
		Ultrasonic system	Air polishing system
A	Standard	Vibration + Water	AirPowder + water
B	No-water mode	Vibration	Air only
C(A+)	Boost mode	5-Level increase in Power	5-Level increase in air pressure
D	Cleaning mode	Water only	Water + Air

Note: Under Boost mode, Power or air pressure would be increased by 5 grades, the maximum is the 12th grade.

After releasing button C, it will automatically return to the previously set grade.


Figure 9.
Schematic diagram of multi-function foot pedal

3.1.2 Standard foot pedal (round)

See description of foot pedal key A in the table above.


Figure 10.
Schematic diagram of multi-function foot pedal

3.2 Ultrasonic system

3.2.1 Scaling

- 1) Turn on the power switch and pick up the ultrasonic handpiece, afterwards the panel will automatically shift to Ultrasonic mode.
- 2) Press the G button to enter Scaling mode.
- 3) Select the tip and tighten it on the handpiece by torque wrench.
- 4) When stepping on the button A, tip vibrates with LED on and water emits; Release the button, vibration and water stop, LED will be off 10 seconds later.
- 5) The handpiece is generally held in a pen-holding position.
- 6) The frequency is extremely high when the machine is on working. It can remove the tartar with the tips softly stroke on the teeth surface like erasing motion. DO NOT stay long or put overmuch pressure on the teeth.
- 7) Vibrating intensity: Adjust the vibrating intensity according to your need. The setting power is level 3, please adjust the vibrating intensity according to the teeth sensitivity and hardness of calculus during clinical application.
- 8) Water volume: Please adjust the volume on the panel.
- 9) During clinical scaling, please keep the side of the tip adapt to the tooth surface horizontally, with light lateral pressure to allow the tip vibrate freely.
- 10) After operation, please keep device working for 30s with water to flush the hand-

piece and tips; Remove the tip and get it sterilized.

3.2.2 Periodontal treatment

- 1) Use a torque wrench to tighten the periodontal treatment tip to the ultrasonic handpiece. Click the P button on panel to enter the periodontal treatment mode.
- 2) The operation and adjustment methods are similar to Scaling mode.

3.2.3 Endodontic irrigation

- 1) Tighten the endo file to the handpiece by endo wrench.
- 2) Click on the E button and enter the Endo mode.
- 3) The default power is level 1 under Endo mode, operator can adjust the power according to the need in clinical treatment.
- 4) Select the appropriate Endo file and slowly place it into the root canal of the patients teeth. Step on the foot pedal to start ultrasonic endodontic irrigation.
- 5) During clinical cleansing, the file should not be pressed too tightly when it is in the root canal
- 6) Do not step on the pedal before the file is in the root canal.
- 7) The suggested power range for endodontic irrigation is level 1-level 5.

3.3 Air polishing system

Use exclusively Xpower powders.

FAST type powders (65µm) must be loaded into FAST powder tank only and used exclusively with PA handpiece.

MILD tupe powders (25/14µm) must be loaded into MILD powder tank only and can be used with both PA and PAP handpieces.

- 1) Add appropriate amount of powder to the powder tank (the amount of powder added should be between "Max" and "Min" marked on the the tank), then tighten the upper cover of the powder tank, and fix the powder tank on the main unit.
- 2) Pick up the Air polishing handpiece, afterwards the panel will automatically shift to Air polishing mode. Select powder type (FAST/MILD) on display.
- 3) Adjust the water volume and air pressure, aim the nozzle at a pool, press the foot button A, and confirm that the nozzle can spray gas, powder and water spray normally before use.
- 4) Before Air polishing treatment, please wear goggles and veil on the face of the patient. Users should wear goggles or a protective mask.
- 5) The handpiece is generally held in a pen-holding position.
- 6) Adjust the water volume and air pressure to an appropriate level. The recommended water volume starts at level 5 and the air pressure starts at level 3. According to the sensitivity of the patient's teeth and hardness of the dental plaque, adjust the water volume and air pressure at any time during the clinical process; increasing the air pressure will enhance cleaning effect, but will weaken the polishing effect; increasing the amount of water will enhance the polishing effect, but will weaken the cleaning effect.
- 7) The nozzle should be aligned with the tooth surface during the cleansing, but avoid direct contact. Keep the nozzle and the tooth surface at a distance of 3-5mm at 30-60 angle. The smaller the angle, the larger the cleaning area; during the cleansing, please Perform a small circular motion on the surface of the tooth.
- 8) The air/powder mixture reflected from the tooth surface should be evacuated by using a strong suction device on the dental unit during treatment.
- 9) After treatment, adjust the water volume to the maximum level, and polish the surface of teeth.

3.4 Cleaning mode

To keep the unit pipeline clean and fresh it is recommended to use Clean&White Water Jet Tablets into water tank.

It is recommended to flush and disinfect the pipeline of the device every day. "Cleaning" mode allows cleaning and disinfection of the pipeline to reduce the accumulation of crystals and the number of bacteria in the pipeline.

- 1) Fill a water tank with distilled water or demineralized water.
- 2) Pick up the ultrasonic handpiece, point the handpiece at the sink, click the Cleaning button on the screen, Step on the pedal to begin cleaning.
- 3) 30s after cleaning, release the pedal to stop.
- 4) After cleaning, put the ultrasonic handpiece back into the bracket. And then, pick up the air polishing handpiece, point the handpiece nozzle at the pool, and click the Cleaning button again, so that the device will automatically blow out the residual powder in the pipeline and release the high pressure gas in powder tank.
- 5) After cleaning for 15 seconds, the device will automatically exit the Cleaning mode. You can also click Purge on the screen to stop cleaning.

3.5 Setting

- 1) Click "ON" / "OFF" to turn on or off the heating function. (According to the actual water output, the water temperature is different).
- 2) Click the "Settings" button at the bottom right corner of the screen to enter the "System Settings" interface
- 3) Click Calibration button, the system will automatically calibrate the air pressure to the value set in the program
- 4) The Bluetooth pedal has been matched to the machine, any subsequent replacement of foot pedal or main unit needs to make new pair.

4 Troubleshooting

4.1 Troubleshooting

Fault	Possible cause	Solutions
The tip does not vibrate and no water comes out after powering on and stepping on the foot switch	Poor contact of power plug	Plug the foot pedal plug well.
	Poor contact of pedal switch	Plug the power supply well
	The fuse is broken	Replace a T1.6AL 250V fuse
	Out of battery under wireless pedal mode	Replace a new battery
	Poor connection under wired pedal mode	Connect the wire and screw it well
The tip does not vibrate and there is water spray after power on and stepping on the foot pedal.	Main unit fault	Contact local distributor or the manufacturer.
	Loose tip	Tighten the tip
	The connection between the tail wire and the circuit board is loose.	Contact local distributor or the manufacturer.
The tip vibrates but there is no water spray after power on and stepping on the foot pedal.	Handpiece fault	Contact local distributor or the manufacturer.
	Cable fault	Contact local distributor or the manufacturer.
After power-off, there is still water spray.	The resistance of the diaphragm pump becomes larger, and the torque of the low water retaining capacity is insufficient	Turn up the water volume.
	Impurity in the solenoid valve	Contact local distributor or the manufacturer.
Handpiece heating	Water line clogging	Contact local distributor or the manufacturer.
	Impurity in the solenoid valve	Contact local distributor or the manufacturer.
Water spray is too small	Water volume is too small	Turn up the water volume
	Device fault	Contact local distributor or the manufacturer.
Weakened tip vibration	Water volume is too small	Turn up the water volume
	Water line clogging	Contact local distributor or the manufacturer.
Outlet water temperature is too high (over 45° C)	Tip is not tightened.	Tighten the tip
	Tip is loose	Tighten the tip
	Broken tip (Note1)	Remove the tip
The Endo file does not vibrate	Thermocouple fault	Contact local distributor or the manufacturer.
There is no air spray and water spray after power on and stepping on the foot pedal.	The nut is not tightened	Tighten the nut
	Loose contact of power supply plug.	Plug the power supply plug well
	Poor connection under wired pedal mode	Connect the wire and screw it well
The nozzle does not spray gas but there is water spray after power on and stepping on the foot pedal.	Out of battery under wireless pedal mode	Replace a new battery
	Nozzle clogging	Dredge the nozzle
	Handpiece clogging	Dredge the handpiece
	Clogging of handpiece tail cord	Remove the tail cord from the main unit, dredge the tail cord or replace it.
There is air flow but no water spray after power on and stepping on the foot pedal.	Solenoid valve fault	Contact local distributor or the manufacturer
	The resistance of the diaphragm pump becomes larger, and the torque of the low water retaining capacity is insufficient	Turn up the water volume
	Impurity in the solenoid valve	Contact local distributor or the manufacturer
Leaking air of powder tank	Water line clogging	Contact local distributor or the manufacturer
	The O-ring on base of powder tank is broken	Unscrew the powder tank base, take out the damaged O-ring, replace the O-ring of the same specification attached, and screw on the powder tank base
	Broken O-ring	Replace the O-ring
	There is powder residue at the thread, so that the screw is not in place.	Remove the residual powder at the thread part.
	The upper cover of powder tank is broken.	Replace the upper cover of powder tank
Water leakage of handpiece	The thread of powder tank is broken so that the screw is not in place	Replace a powder tank
	Broken O-ring of handpiece	Replace the O-ring
The air powder scaling efficiency is reduced.	The powder in tank is not enough	Add powder to the tank
	Powder residue in pipe, handpiece, or nozzle passage	Clean the passage with a fine needle and blow it off with compressed air
Touch panel pops up prompt message	No powder tank	Check the powder tank and reinstall
	Air pressure is too high or too low	Increase the external air pressure input, Input air pressure range requirements: 5.5bar-7.5bar(0.55MPa-0.75MPa)
	Pick up two handpieces at a time	Before stepping on the pedals, please pick up only one handpiece
	No handpiece has been picked up	Before stepping on the pedals, please pick up only one handpiece
	The handpiece does not match the mode configuration	Please select the correct handpiece and mode

Note: if the problems cannot be solved, please contact local distributor or the manufacturer.

4.2 Notes

Note1 If the tip is surely tightened and there is water mist spray, the tip is considered to be damaged with the following phenomena:

- 1) The vibration intensity of the tip and the degree of water atomization are significantly weakened.
- 2) The tip makes a harsh click sound during operation.

5 Cleaning, disinfection, and sterilization

5.1 Please conform to the recommendations of the attachment "Reprocessing Instructions of Cleaning, Disinfecting and Sterilizing" delivered with your product regarding the procedure of cleaning, disinfecting, sterilizing and packing of the components. Follow present-day regulations enforced in the country about reprocessing. The use of scouring powder or an abrasive sponge will damage its surface.

5.2 The handpiece, tips, torque wrench, endo wrench, LED lamp and Light conductor (the handpiece with LED lamp) can be sterilized.

- Recommended sterilization cycles of handpiece: 300 cycles
- Recommended sterilization cycles of torque wrench and endo wrench: 300 cycles

Notice:

- a) Clean the handpiece with compressed air before sterilization.
- b) Be sure that the tip has been unscrewed from the handpiece and it cannot be sterilized with others.
- c) Please notice whether the outer of the handpiece is damaged during the treatment and sterilization. Don't smear any protective oil on the surface of handpiece.
- d) There are some waterproof "O" rings at the end of the handpiece. Please lubricate them with dental lube frequently, as sterilization and repeated pulling and inserting will reduce their using life. Change a new one once it is damaged or worn excessively.

e) The following sterilizing methods are forbidden:

- Boil in water.
- Dip in iodine, alcohol and glutaraldehyde.
- Bake in oven or microwave oven.

Notice: We are not responsible for any damage caused in the above items.

6 Transport, storage and maintenance

6.1 Transport

- 6.1.1 Excessive impact and shake should be prevented during transport. Lay it care fully and lightly.
- 6.1.2 Do not put it together with dangerous goods during transport.
- 6.1.3 Avoid being exposed to sun, rain, and snow during transport.

6.2 Storage

- 6.2.1 The device should be handled carefully and lightly. Be sure that it is far from the vibration, and installed or kept in a cool, dry, and ventilated place.
- 6.2.2 Do not store the machine together with articles that is poisonous, combustible, caustic, or explosive.
- 6.2.3 This machine should be stored in a room where the relative humidity is 10% - 93%, atmospheric pressure is 70kPa - 106kPa, and the temperature is -20°C - +5°C.
- 6.2.4 When the device is not in use, turn off the power supply and unplug the power plug. If it is not used for a long time, it should be energized and connect to water and air once a month for five minutes.

6.3 maintenance

6.3.1 Air filter

- 1) When water accumulate in the filter, turn the knob at the bottom of the filter counterclockwise to drain the water, and tighten the knob clockwise.
- 2) Replacement of the filter element: Use a filter wrench to unscrew the transparent cover of the air filter, then use the wrench to unscrew the black nut at the lower end of the filter element, remove the white filter element and discard it into the trash can, re- place it with a new filter element, and reinstall the black nut and transparent shell. It is recommended to replace the filter element every 24 months, and the spare filter element is included in the accessory.


Figure 11. Schematic diagram of filter element replacement

7 Environmental Protection

The instrument contains no hazardous components and can be disposed of and destroyed in accordance with the relevant local regulations.

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○ : Indicates that the content of the toxic substance in all homogeneous materials of the part is below the limit requirement stipulated in SJ/T-11363-2006 Limit Requirements for Toxic and Hazardous Substances in Electronic Information Products.

X:Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceed the limit requirement specified in SJ/T-11363-2006. Please dispose according to the local laws or consult with dealer from whom you purchased it about waste disposal.

8 After service

- 8.1** After the equipment is sold, the manufacturer will be responsible for quality problem according to the warranty card. For specific items, please refer to the warranty instructions in the warranty card.
- 8.2** This product does not contain self-maintaining parts. All maintenance, adjustment, calibration, and modification of technical parameters of the product can only be carried out by the technicians or special repair shops. If the customer needs to repair by himself, the manufacturer can provide circuit diagrams, component lists, legends, calibration rules, or other information necessary to help the user's qualified technicians repair the equipment parts designated by the manufacturer, but the manufacturer will not bear the consequences arising therefrom.
- 8.3** The user must use the original accessories, please contact your local dealer or the manufacturer to purchase.
It is forbidden to use accessories of other brands to avoid damage to the equipment or other dangers.
- 8.4** After the handpiece, tips and other accessories are damaged, users should not repair them by themselves.
Please purchase new parts and replace them before use. If you need relevant information, please contact the manufacturer.
- 8.5** Repair service only valid for the EU territory, for countries outside the EU please contact your dealer. All the equipment distributed by us that will need technical interventions during or outside the warranty period, must be sent directly to the company Vi.Vi.Med srl who will manage the repair and invoice.
- VI.VI.MED Srl - Technical Service**
Via Tevere 11 - 51018 Pieve a Nievole (PT) - Italy
Tel. +39 0572 520910 - fax. +39 0572 520807 - Mobile: +39 347 5250335
e-mail: info@vivimed srl - web www.vivimed srl.com

9 EMC-Declaration of conformity

9.1 Instructions for use

The ME EQUIPMENT or ME SYSTEM is suitable for hospitals and dental clinics.

Warning: Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Note: The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

List of all cables

NO	Name	Length	Shielded or not	Detachable or not	Note
1	Power cord	2.0m	No	Yes	/
2	Foot pedal cord	2.8m	No	Yes	/
3	Ultrasonic handpiece cord	2.0m	No	No	/
4	Air polishing handpiece cord	2.0m	No	No	/

Replaceable accessories please refer to 1.6 replacement instruction for main accessories.

Performance of the me equipment

Xpower Ultrasonic Scaler & Air Polisher realizes the ultrasonic cleaning function by vibration of ultrasonic handpiece. According to different applications, "G", "P", "E" modes can be selected, "A" mode Air polishing function is adjusted by controlling the air pressure to achieve adjustment of air polishing power; "C" mode is the cleaning mode of the ultrasonic handpiece or air polishing handpiece. After the operation, the ultrasonic handpiece pipeline and the air polishing handpiece pipeline need to be cleaned. When the me equipment essential performance is lost or degraded due to em disturbances, the doctor should immediately stop using it to ensure that there is no treatment error. And then remove the source of disturbances or adjust the direction or position of me equipment to ensure me equipment can be used in normal performance condition.

9.2 Technical description

9.2.1 All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the excepted service life.

(1) Portable and mobile RF communications equipment may affect the performance of equipment, use of equipment should be avoided strong electromagnetic interference, and do not closer to mobile phone, microwave oven, etc.

(2) Use of equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

(3) Except for the cables sold by manufacturers of as spare parts of internal components, the use of accessories and cables other than those specified or provided by the manufacturer may result in increased electromagnetic emissions or decreased

electromagnetic immunity of this equipment and result in improper operation.

(4) Use of accessories, transducers and cables other than those specified or provided by the manufacturer together with equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

9.2.2 Guidance and manufacturer's declaration -electromagnetic emissions and Immunity

Table 1

Guidance and manufacturer's declaration - electromagnetic emissions	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class A
Harmonic emissions IEC 61000-3-2	Not Applicable
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable

Table 2

Guidance and manufacturer's declaration - electromagnetic Immunity		
Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±2 kV power supply lines Not applicable 100 kHz repetition frequency
Surge IEC 61000-4-5	±0.5 kV, ±1 kV differential mode ±0.5 kV, ±1 kV, ±2 kV common mode	±0.5 kV, ±1 kV differential mode Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle	0 % UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz/30 A/m	3 V 0,15 MHz – 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz

NOTE UT is the a.c. mains voltage prior to application of the test level.

Table 3

Guidance and manufacturer's declaration - electromagnetic Immunity						
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band (MHz)	Service	Modulation	IEC 60601-1-2 Test Level (V/m)	Compliance level (V/m)
		385	380–390	TETRA 400	Pulse modulation 18 Hz	27
450		430–470	GMRS 460,FRS 460	FM ± 5 kHz deviation 1 kHz sine	28	28
710		704–787	LTE Band 13,17	Pulse modulation 217 Hz	9	9
745						
780						
810		800-900	GSM 800/900,TETRA 800, IDEN 820,CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28	28
870						
930						
1720						
1845		1700-1990	GSM 1800;CDMA 1900; GSM 1900;DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	28	28
1970						
2450		2400-2750	Bluetooth, WLAN,802.11 b/g/n, RFID 2450,LTE Band 7	Pulse modulation 217 Hz	28	28
5240	5100–5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9	9	
5500						
5785						

Table 4

Guidance and manufacturer's declaration - electromagnetic Immunity				
Radiated RF IEC61000-4-39 (Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields)	Test Frequency	Modulation	IEC 60601-1-2 Test Level (A/m)	Compliance level (V/m)
	30 kHz	CW	8	8
	134,2 kHz	Pulse modulation 2.1 kHz	65	65
	13,56 MHz	Pulse modulation 50 kHz	7,5	7,5

Attachment Reprocessing instructions of cleaning, disinfecting and sterilizing

1. Beginning work

1.1 Please read these operating instructions carefully as they explain all the most important details and procedures. Please pay special attention to the safety precautions. Always keep this instruction close at hand.

1.2 To prevent injury to people and damage to property, please heed the corresponding directives.

1.3 The instructions in this manual are only applicable to the product which it was delivered with.

2. Introduction

- 2.1 These reprocessing instructions provide instructions for cleaning, disinfection, sterilization and packaging of manufacturer reusable products intended to be reprocessed in medical facilities.
- 2.2 The goal of reprocessing reusable products is to reduce bioburden and to achieve sterility of those products in order to eliminate the risk of product reuse related infection. Decisions regarding cleaning, disinfecting or sterilizing manufacturer's medical and dental instruments are based on the potential risk of infection associated with their use.
- 2.3 It is recommended to use steam sterilization.
- 2.4 Remember that sterilization or high-level disinfection cannot be achieved unless the elements of the assembly are cleaned first.
- 2.5 If you find that the reprocessing instructions from the manufacturer seem to be inadequate, please inform manufacturer about those inadequacies.
- 2.6 We encourage you to report adverse events related to device reprocessing. Report such events directly to manufacturer.

3. Reprocessing - instructions for reusable products

- 3.1 The instructions are binding for the reprocessing of all reusable products (Here after called "products") of manufacturer. When necessary, additional product-specific instructions are included with the product to provide additional information.
Important: Before use, carefully read the operating instructions of the manufacturer instrument and devices with which the product will be used.
- 3.2 Reusable products must be cleaned, disinfected and sterilized prior to first use. Reprocessing procedures have only limited implications to this device. The limitation of the numbers of reprocessing procedures is therefore determined by the function / wear of the device. From the processing side there is no maximum number of allowable reprocessing. The device should no longer be reused in case of signs of material degradation. In case of damage the product should be reprocessed before sending back to the manufacturer for repair.

4. Preparation - basic principles

- 4.1 It is only possible to carry out effective sterilization after the completion of effective cleaning and disinfection. Please ensure that, as part of your responsibility for the sterility of products during use, only sufficiently validated equipment and product-specific procedures are used for cleaning/disinfection and sterilization, and that the validated parameters are adhered to during every cycle.
- 4.2 Please also observe the applicable legal requirements in your country as well as the hygiene regulations of the hospital or clinic. This applies especially with regard to the additional requirements for the inactivation of prions.

5. Preparation at the point of use

Disconnect product. Remove gross soiling of the products with cold water (<40° C) immediately after use. Don't use a fixating detergent or hot water (>40° C) as this can cause the fixation of residuals which may influence the result of the reprocessing process. Store the products in a humid surrounding.

6. Transportation

Safe storage and transportation to the reprocessing area to avoid any damage and contamination to the environment.

7. Preparation for decontamination

The products must be reprocessed in a disassembled state, as far as possible.

8. Pre-cleaning

Do a manual pre-cleaning, until the products are visually clean. Submerge the products in a cleaning solution and flush the lumens with a water jet pistol with cold tap water for at least 10 seconds. Clean the surfaces with a soft bristle brush.

9. Cleaning

Regarding cleaning/disinfection, rinsing and drying, it is to distinguish between manual and automated reprocessing methods. Preference is to be given to automated reprocessing methods, especially due to the better standardizing potential and industrial safety.

Automated Cleaning:

- Use a washer-disinfector (WD) meeting the requirements of the ISO 15883 series.
Put the products into the machine on a tray. Connect the products with the WD by using suitable adapter and start the program:
4 min pre-washing with cold water (<40° C)
Emptying
5 min washing with a mild alkaline cleaner at 55° C
Emptying
3 min neutralising with warm water (>40° C);
Emptying
5 min intermediate rinsing with warm water (>40° C)
Emptying

The automated cleaning processes have been validated by using 0.5% neodisher MediClean forte (Dr. Weigert).

⚠ Acc. to EN ISO 17664 no manual reprocessing methods are required for these devices. If a manual reprocessing method has to be used, please validate it prior to use.

10. Disinfection

Automated Thermal Disinfection in washer/disinfector under consideration of national requirements in regards to A0 value (see EN 15883).
A disinfection cycle of 5 min disinfection at 93° C has been validated for the product to achieve an A0 value of 3000.

11. Drying

Automated Drying:

Drying of outside of products at 40°C , 5 min through drying cycle of washer/disinfector. If needed, additional manual drying can be performed through lint free towel. Insufflate cavities of products byusing sterile compressed air.

12. Functional testing, maintenance

Visual inspection for cleanliness of the products and reassembling if required. Functional testing according to instructions of use. If necessary, perform reprocessing process again until products is visibly clean.

Before packaging and autoclaving, make sure that the products have been maintained acc. to manufacturer's instruction.

13. Packaging

Pack the products in an appropriate packaging material for sterilization. The packaging material and system refer to EN ISO 11607.

14. Sterilization

Sterilization of products by applying a fractionated pre-vacuum steam sterilization process (according to EN 285/EN 13060/EN ISO 17665) under consideration of the respective country requirements.

Minimum requirements: 3 min at 134° C (in EU: 5 min at 134° C)

Maximum sterilization temperature: 138° C

Drying time:

For steam sterilization, we recommend a drying time of 15 to 40 minutes. Choose a suitable drying time, depending on the autoclave and load. Refer to the autoclave's instructions for use.

After sterilization:

- Remove the product from the autoclave.
- Let the product cool down at room temperature for at least 30 minutes. Do not use additional cooling.

Check that the sterilization wraps or pouches are not damaged.

⚠ Flash sterilization is not allowed on lumen products.

⚠ The manufacturer assumes no responsibility for the use of other sterilization procedures (e.g. ethylene oxide, formaldehyde and low temperature plasma sterilization).

In such cases,

please observe the respective valid standards (EN ISO 14937/ANSI AAMI ISO 14937 or the procedure specific standard) and verify the suitability and effectiveness in principle of the procedure (if necessary, including investigations on sterilizing agent residue), taking into account the specific product geometry as part of the validation.

- Maximum sterilization temperature 138° C

15. Storage

Storage of sterilized products in a dry, clean and dust free environment with a relative humidity of 10% to 93%, an atmospheric pressure of 70KPa to 106KPa, and a temperature of -20° C to +55° C; refer to label and instructions for use.

After sterilization, the product should be packaged in a medical sterilization bag or a clean sealing container, and stored in a special storage cabinet. The storage time should not exceed 7 days. If it is exceeded, it should be reprocessed before use.

16. Service life

The products have been designed for a large number of sterilization cycles. The materials used in their manufacture were selected accordingly. However with every renewed preparation for use, thermal and chemical stresses will result in aging of the devices. If the number of permissible re-sterilization cycles is restricted, this will be pointed out in the product specific instructions.

⚠ The use of ultrasound baths and strong cleaning and disinfection fluids (alkalizing pH>9 or acid pH<5) can reduce the life span of devices. The manufacturer accepts no liability in such cases.

⚠ The devices may not be exposed to temperatures above 138° C.

It is the duty of the user to ensure that the reprocessing processes including resources, materials and personnel are capable to reach the required results. State of the art and often national law requiring these processes and included resources to be validated and maintained properly.

10 Disposal and Scrapping

The product doesn't contain battery or toxic substances. And there are no components which should be removed specially from the main unit for disposal and scrapping.

After the device is out of its service life, you must not discard it in domestic household waste. Please comply with the Waste Electrical and Electronic Equipment (WEEE) directives and the medical waste disposal regulations of your country.

Torque wrench, endo wrench and handpiece, those which could easily contact to the biological sources and cause biological hazards, shall be detached from the main unit and reprocessed according to the reprocessing treatment in Chapter 5 before the disposal and scrapping.

11 Statement

Refine reserves the right to modify the product technology, accessories, instruction manual and product packaging content at any time without notice. The product is subject to the actual product, and the pictures are for reference only. The final interpretation right belongs to Guilin Refine Medical Instrument Co., Ltd.

Shelf life: 10 years, the date of manufacture see product label.

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MedNet EC-REP C IIb GmbH
MediMap Ltd