OPTICORE DUAL BOND

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GB OPTICORE DUAL BOND ADHESIVE

Recommendation for Use Light cure single component adhesive

OPTICORE DUAL BOND ADHESIVE is a simple to use one bottle light cure single component priming and dentin/enamel adhesive in one. It is designed for strong bonding of composites, compomers and metals to enamel and dentin and nonprecious and precious metals. OPTICORE DUAL BOND ADHESIVE can be also used for priming the root canals before filling of the root canals or cementation of endodontic posts.

OPTICORE DUAL BOND ADHESIVE is further designed as a primer for indirect adhesive restorations, e.g. ceramic, metal and composite inlays, onlays, veneers, crowns and bridges with a self or dual curing composite cement. Its strong adhesion to etched enamel or dentin proceeds on principles similar to that occurring with glassionomer cements. Good, long lasting adhesive strength and good biocompatibility are attained by methacrylate grafted polycarboxylic acids.

OPTICORE DUAL BOND ADHESIVE is compatible with all current brands of visible light cure composite restorative materials. It is ethanol-based and hydrophilic. Therefore OPTICORE DUAL BOND ADHESIVE can be used on slightly moist dentin surfaces (wet bonding technique). In the rare cases a dual or self cure priming and bonding system is recommended OPTICORE DUAL BOND ADHESIVE can get dual cure which is applied as a 1:1-mixture of OPTICORE DUAL BOND ADHESIVE and OPTICORE DUAL BOND ACTIVATOR. The activator is available separately. For these applications refer to recommendations for use of OPTICORE DUAL BOND ACTIVATOR.

1. Priming of root canals

1.1. Preparation

Use of a rubber dam to isolate the tooth is strongly recommended.

Prepare and clean the root canal with hydrogen peroxide (Opticore Canal Cleaner) or sodium hypochlorite solution.

For fixing of endodontic posts with compo-

site cements: Following the directions of the selected post manufacturer, prepare the canal to accept the appropriate length post. Trial fit the post and adjust for length, if necessary. Remove the post and prepare the canal for the endodontic post cementation.

1.2. Priming of root canals

Because the surfaces of the root canal surfaces are not accessible by a brush, use a syringe for the application of OPTICORE DUAL BOND ADHESIVE and flush it into the prepared root canals. The flushing must go on over a time of 60 seconds discontinuous in 4 cycles. All the dentin surfaces must kept wet with the primer over the specified time by generously brushing the accessible surfaces with a brush.

An excess of the primer in the canals is to be removed with soft tissue paper. Dry the root canal cautiously with oilfree air for 15 seconds and cure the OPTICORE DUAL BOND ADHESIVE coating with a dental halogen light unit for 20 seconds.

1.3. Cement placement

Refer to manufacturer instructions for placement and curing of the composite cements.

2. Bonding of Light Cure Composite and Composer Filling Restorations

2.1. Isolation

Rubber dam is the recommended method of isolation.

2.2. Cavity Preparation

Clean the tooth with flour of pumice and water prior to preparation. Prepare the cavity with minimal tooth reduction. Margins should have a slight (0.5 - 1.0 mm) bevel placed in the enamel to increase the surface area for greater bond strength.

2.3. Pulp Protection

Cavity floor of deep excavations should be covered with a thin layer of hard setting calciumhydroxide material.

2.4. Enamel and Dentine Conditioning

Recommended is the total etch technique. Apply OPTIMA BLUETCH onto the enamel and dentine surfaces beginning with the enamel bevels. Leave the etching gel in place for 20 seconds. On primary teeth and teeth high in fluoride, a 60 second etch is recommended. Rinse for 20 seconds with water. Dry it in a water and oil-free airstream, but do not desiccate. A slightly wet dentine surface is important for the function of OPTICORE DUAL BOND ADHESIVE. The etched enamel bevel should have a chalky white appearance.

Etching precaution:

It is essential, that etched areas are not contaminated by anything. If contamination occurs, re-etch, rinse with water and dry as above. Avoid etching gel contact with oral soft tissues, eyes and skin. If accidental contact occurs, flush immediately with copious amounts of water.

For the typical indications of compomers etching is mostly not recommended, but it improves considerably the bonding with tooth material.

2.5. Application of OPTICORE DUAL BOND ADHESIVE

Apply OPTICORE DUAL BOND ADHESIVE generously with a brush onto the enamel and dentin surfaces for 30 seconds with agitation. Dry cautiously with oil free air for 15 seconds. Do not desiccate the dentin. Cure the OPTICORE DUAL BOND ADHESIVE coating by exposing its entire area to a dental halogen light unit for 20 seconds before application of a second layer of OPTICORE DUAL BOND ADHESIVE.

Notes: Do not rinse off the OPTICORE DUAL BOND ADHESIVE! If not used immediately, place dispensed OPTICORE DUAL BOND ADHESIVE in subdued light to prevent premature polymerization by incident light. The OPTICORE DUAL BOND ADHESIVE will not

selfcure.

2.6 Application of a Second Layer of OPTICORE DUAL BOND ADHESIVE

Apply again the OPTICORE DUAL BOND ADHESIVE generously with a brush onto the adhesive surfaces for 30 seconds with agitation.

Thoroughly dry again the primed areas with oil free air for

15 seconds and cure the OPTICORE DUAL BOND ADHESIVE coating by exposing its entire area to a dental halogen light unit for 20 seconds before placement of a light cure composite.

2.7. Restorative Placement, Cure and Finishing

Refer to manufacturers instructions for placement, curing and finishing of light cure restorative materials.

3. Sealing of the Cavity before Amalgam Filling Restorations

3.1. Application of Light Cure OPTICORE

DUAL BOND ADHESIVE to Dentin/Enamel After isolation of tooth, preparation of cavity, and in case of need pulp protection according point 1.1., 1.2. and 1.3., apply OPTICORE DUAL BOND ADHESIVE in two layers by the same way as in 1.5 and 1.6 described.

It is essential that the primed areas are dry and contaminant free for the amalgam application. Because light curing under amalgam restorations is impossible, and the primer does not selfcure it is essential to cure the coating by exposing its entire area to a dental halogen light unit for

20 seconds before application of the amalgam filling. Alternatively, please use the OPTICORE DUAL BOND ADHESIVE in the dual cure mode (in addition with OPTICORE DUAL BOND ACTIVATOR).

3.2. Application and Condensing of Amalgam Apply amalgam according instructions of the manufacturer.

4. Application of OPTICORE DUAL BOND ADHESIVE on Precious and Nonprecious Metals and Ceramics for repairing of crowns and bridges

4.1. Treatment of the surfaces of the prosthetic works

Rough surfaces of metallic or fractured ceramic surfaces with a diamond bur (40µm) or with suitable intraoral sand blasting. Treat the roughen surfaces with OPTIMA BLUETCH, then rinse with plenty of water and dry afterwards with oil and waterfree air. Prevent contact of the etchant with skin eyes and gingival tissue.

Optional for ceramics:

If the prosthetic work can be done extraorally, ceramics can be etched with 1% hydrofluoric acid, but caution: hydrofluoric acid is strong corrosive. Do never use it intraorally.

It is recommended to prime the prosthetic work first with a suitable silane primer prior to application of OPTICORE DUAL BOND ADHESIVE. Refer to instruction for use of the silane primer.

4.2. Application of OPTICORE DUAL BOND ADHESIVE

Apply OPTICORE DUAL BOND ADHESIVE generously with a brush onto the metallic surfaces and dry these surfaces with air that is oil and water-free for 15 seconds. Then cure the primer with a dental halogen light unit for 20 seconds and start the repair work.

5. Reparations of fillings with amalgam

5.1. Treatment of the amalgam surfaces

Rough metallic surfaces with a diamond bur (40µm) or with suitable intraoral sand blasting. Areas of enamel & dentine are to etch by the same way as in 1.4. described

5.2. Application of OPTICORE DUAL BOND ADHESIVE

Apply OPTICORE DUAL BOND ADHESIVE in two layers according Chapter 1.5. and 1.6.

5.3. Application of a composite

For repairing the amalgam filling, apply a light curing composite according manufacturer instructions.

6. Storage of OPTICORE DUAL BOND ADHESIVE

Do not store above 25 °C (78 °F)! Avoid storage in direct sunlight. Do not use after expiration date.

7. Technical Data

Time for light cure with a dental halogen unit: 20 sec

Setting-time of the 1:1-mix of OPTICORE DUAL BOND ADHESIVE / OPTICORE DUAL BOND ACTIVATOR for autocure (without any light; air excluded): ca. 3 min (refer to recommendation for use for OPTICORE DUAL BOND ACTIVATOR).

OPTICORE DUAL BOND ACTIVATOR Recommendations for Use

OPTICORE DUAL BOND ACTIVATOR is designed as a product for additional and/or extended applications of OPTICORE DUAL BOND ADHESIVE for all the special applications, where a dual or self cure primer and bonding system is strongly recommended.

1. Adhesive cementing of indirect restorations

1.1. Isolation

Rubber dam is the recommended method of isolation.

1.2. Cavity Preparation

Remove any debris and stain prior to preparation. Prepare the cavity with minimal tooth reduction. Margins should have a slight (0.5 - 1.0 mm) bevel placed in the enamel to increase the surface area for greater bond strength.

1.3. Pulp Protection

Cavity floor of deep excavations should be covered with a thin layer of hard setting calcium-hydroxide material.

1.4. Enamel and Dentine Conditioning

Recommended is the total etch technique. Apply OPTIMA BLUETCH onto the enamel and dentine surfaces beginning with the enamel bevels. Leave the etching gel in place for 20 seconds. On primary teeth and teeth high in fluoride, a 60 second etch is recommended.. Rinse for 20 seconds with water. Dry it in a water and oil-free airstream, but do not desiccate. A slightly wet dentine surface is important for the function of OPTICORE DUAL BOND ADHESIVE. The etched enamel bevel should have a chalky white appearance.

Etching precaution:

It is essential, that etched areas are not contaminated by anything. If contamination occurs, re-etch, rinse with water and dry as above. Avoid etching gel contact with oral soft tissues, eyes and skin. If accidental contact occurs, flush immediately with copious amounts of water.

Preparation and application OPTICORE DUAL BOND ADHESIVE / OPTICORE DUAL BOND ACTIVATOR-mix

After cavity preparation, dispense 1 drop of OPTICORE DUAL BOND ADHESIVE on the mixing pad. Add 1 drop of OPTICORE DUAL BOND ACTIVATOR, and mix all together for 15 seconds under subdued light.

1.6. Application to dentin/enamel

Apply immediately the OPTICORE DUAL BOND ADHESIVE / OPTICORE DUAL BOND ACTIVATOR-mix generously with a brush onto the whole area of the enamel and dentin surfaces for 30 seconds with agitation and apply an uniform coating on the whole area of dentin and enamel. Re-apply the primer mixture as needed to insure that the dentin surfaces are kept wet with the primer mixture for the specified time. Dry cautiously the primed areas with oil free air for 15 seconds. The entire area of the coating may be exposed to a suitable light curing unit for 20 seconds before application of a second layer according to 1.7.

1.7. Application of a second layer

Apply again the dual curing 1:1-mixture of OPTICORE DUAL BOND ADHESIVE and OPTICORE DUAL BOND ACTIVATOR generously with a brush onto the adhesive surfaces for 30 seconds with agitation. Dry cautiously again the primed areas with oil free air for 15 seconds and - if desired - expose it to a light curing unit for 20 seconds. It is essential that the primed areas are dry and contaminant free for the application of the indirect restorations. The dual cure system will autocure within ca. 3 minutes after application under the indirect restoration. The indirect restoration is cemented with a dual- or self cure composite cement.

1.8. Application of Indirect Restorations

Apply the indirect restorations according to the instructions of the manufacturer.

2. Bonding of Self Cure or Dual Cure Composites and Compomers

2.1. Preparation and Application of OPTICORE DUAL BOND ADHESIVE / OPTICORE DUAL BOND ACTIVATOR-mix

After cavity preparation according 1.1., 1.2, 1.3, and 1.4., dispense 1 drop of OPTICORE DUAL BOND ADHESIVE on the mixing pad. Add 1 drop of OPTICORE DUAL BOND ACTIVATOR, and mix all together for 15 seconds under subdued light.

2.2. Application to dentin/enamel

Apply this mixture according Chapter 1.6. The coating of OPTICORE DUAL BOND ADHESIVE / OPTICORE DUAL BOND ACTIVATOR may be exposed to a suitable light curing unit for 20 seconds prior to application of a second layer according Point 2.3.

2.3. Application of a second layer

Apply again the dual cure 1:1-mixture of OPTICORE DUAL BOND ADHESIVE and OPTICORE DUAL BOND ACTIVATOR generously with a brush onto the adhesive surfaces for 30 seconds with agitation.

Dry cautiously again the primed areas with oil free air for 15 seconds and - if desired - expose it to a light curing unit for 20 seconds.. It is essential that the primed areas are dry and free of contamination for the application of the self cure or dual cure composites and compomers.

The dual cure system will autocure from time of application and meantime restorative placement within ca. 3 minutes.

2.4. Restorative placement, cure and finishing Refer to manufacturers instructions for placement, curing and finishing of dual cure or self cure restorative materials.

3. Storage

Do not store above 25 °C (78 °F). Do not use after expiration date.