## **I** OPTICORE ETCH BOND Adesivo duale auto mordenzante

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**GB** **OPTICORE ETCH BOND**  
**Self-etching dual cure adhesive**

Rev0 - Dec 22, 2008

OPTICORE ETCH BOND is a simple to use **self-etching dual cure adhesive**. It is designed for strong bonding of **composites, compomers and resin modified glass ionomer cements to enamel and dentin**.

OPTICORE ETCH BOND consists of the components PART A and PART B that were mixed before application. It is bonding to dual cure, self-cure and light cure composite restorative materials.

**APPLICATION**

**1. Isolation**

Rubber dam is the recommended method of isolation.

**2. Cavity preparation**

Clean the tooth with flour of pumice and water prior to preparation. Prepare the cavity with minimal tooth reduction. Margins should have a slight **(0.5 - 1.0 mm)** bevel placed in the enamel to increase the surface area for greater bond strength.

**3. Pulp protection**

Cavity floor of deep excavations should be covered with a

thin layer of hard setting calciumhydroxide material.

**4. Application of OPTICORE ETCH BOND**

**ETCH BOND STANDARD KIT** - One drop of OPTICORE ETCH BOND PART A and one drop of OPTICORE ETCH BOND PART B were combined in a mixing pallet and mixed for 5-10 seconds.

**ETCH BOND SINGLE MIX EASY** - Pierce the supplied microbrush **vertical** through the aluminium foils of the SingleMix Easy unit (up to the ruffles of the microbrush) and mix the two components for 10 seconds.

Apply the homogeneous mixture generously with a brush onto the slightly wet enamel- and dentin surfaces for 30 seconds with agitation. The material should build a homogeneous layer. Air thin for 10 seconds to remove the volatile components and to disperse the adhesive. Then light cure with a dental halogen light unit for 20 seconds and place the restorative material. For a maximal adhesion it is **strongly recommended to use the light cure mode**.

If light cure is absolutely impossible, the adhesive will also cure in the autocure mode. After application of the

adhesive air thin to **remove all volatile components**. There must remain a sticky layer. Then apply the restorative material. The setting time of the restoration is 3-5 minutes in the autocure mode.

Notes: *Do not interchange lids of the bottles, because this can lead to a cross-contamination of the liquids.*

**5. Restorative Placement, Cure and Finishing**

Refer to manufacturer instructions for placement, curing and finishing of the used restorative materials.

**STORAGE**

Do not store above 25 °C (78 °F)! Avoid storage in direct sunlight. Do not use after expiration date.

**ATTENTION**

This product is intended for professional dental use. Selling of this product outside UE could be subject to limitations.