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## GB OPTICORE Z Dual Curing Flowable Hybrid Composite for Core Build-ups and Cementation of Posts

OPTICORE Z is a **dual cure** high radiopaque flowable **hybride composite** material with zirconium-dioxide particles for **core-build-up** and **cementation of posts**. Due to its excellent mechanical properties final crown preparation can be carried out more precise. Delivered AUTOMIX-cartridges, OPTICORE Z can be easily dispensed and applied directly. OPTICORE Z exhibits a short setting time without high heat generation. The dual cure properties enables the dentist, also to carry out cementations and core build-ups in cases where a light cure cannot be guaranteed to be sufficient.

### Preparing the cartridge

Remove the cap of the cartridge and throw it away (**do not use it again!**). It is replaced by a special 1:1 mixing cannula. Turn the cannula 90° until it locks in position. Fix the oral tip to the mixing cannula and the ENDO TIP onto the oral tip. The Material is now ready for application.

*Nota: Discard the first 2-3 mm of the extruded material. This must be done for each new mix.*

Store used cartridge with fixed used mixing cannula in the dark.

**The working time (23°C (74°F)) of OPTICORE Z in the self-cure mode is 1:30 minutes from start of mixing.**

## 1 POST CEMENTATION

### 1.1 Isolation

Use of a rubber dam to isolate the tooth is strongly recommended.

### 1.2 Root Canal Preparation

Refer to directions of the selected post manufacturer. (see diagram 1, 2, 3)

Remove all existing old restorations and decay from the tooth. Prepare and clean the root canal with e.g. sodium hypochlorite solution, rinse and remove excess solution from the canal with a soft paper tip. (see diagram 4)

Condition and apply a bonding agent according to the manufacturer instructions. The total etch technique with OPTIETCH and subsequent application of a dual cure bonding agent (e.g. OPTICORE DUAL BOND and ACTIVATOR) is recommended. (see diagram 5a, 5b)

Otherwise use a self etching bonding agent (see diagram 5):

One drop of OPTICORE ETCH BOND PART A and one drop of OPTICORE ETCH BOND PART B were combined in a mixing pallet and mixed for 5-10 seconds.

*Nota: Do not interchange lids of the bottles, because this can lead to a cross-contamination of the liquids.*

Apply the homogeneous mixture generously with a brush onto the slightly wet root canal walls for 15 seconds with agitation. Repeat procedure 1 - 2 times. All the dentine surfaces must kept wet with the primer over the specified time.

Remove excess OPTICORE ETCH BOND with dry paper points. Dry the root canal cautiously with oilfree

air for 15 seconds to remove all volatile components and to disperse the adhesive to an even layer. If not using translucent OPTICORE POST, light cure all areas that are available for a dental curing unit for 20 seconds.

### 1.3 Post cementation

Prepare the selected post according to manufacturer directions.

Using the ENDO TIP, OPTICORE Z is applied into the prepared root canal and onto the post. Seat the post careful into the canal and maintain firm pressure until the post is seated. OPTICORE Z self-cures within **3:30 minutes**. For post stabilization light cure the coronal part of the cemented post for **20 seconds**. (see diagram 6, 7)

Remove the ENDO TIP from the oral tip and proceed with the core-build-up procedure.

## 2 CORE BUILD-UP

### 2.1 Isolation

Use of a rubber dam to isolate the tooth is strongly recommended.

### 2.2 Cavity preparation

Remove all existing old restorations and decay from the tooth. If necessary place any pins or posts. Refer to directions of the selected post manufacturer.

### 2.3 Pulp protection

Cavity floor of deep excavations should be covered with a thin layer of hard setting calcium-hydroxide material.

### 2.4 Application of OPTICORE ETCH BOND

Condition and apply a bonding agent according to the manufacturer instructions. The total etch technique with OPTIETCH and subsequent application of a dual cure bonding agent (e.g. OPTICORE DUAL BOND and ACTIVATOR) is recommended.

Otherwise use a self etching bonding agent:

One drop of OPTICORE ETCH BOND PART A and one drop of OPTICORE ETCH BOND PART B were combined in a mixing pallet and mixed for 5-10 seconds.

*Nota: Do not interchange lids of the bottles, because this can lead to a cross-contamination of the liquids.*

Apply the homogeneous mixture generously with a brush onto the slightly wet enamel and dentin surfaces for 30 seconds with agitation. The material should build a homogeneous layer. Air thin for 10 seconds to remove the volatile components and to disperse the adhesive. Then light cure with a dental halogen light unit for 20 seconds and place the restorative material. For a maximal adhesion it is **strongly recommended to use the light cure mode.**

If light cure is absolutely impossible, the adhesive will also cure in the autocure mode. After application of

the adhesive air thin to **remove all volatile components**. There must remain a sticky layer.

*Nota: It is essential that the primed dentine and enamel surfaces are dry and contaminant free for the application of OPTICORE Z.*

### 2.5 Application

Place the mixing cannula directly into the preparation and press out the paste.

OPTICORE Z is automatically mixed when dispensed with slight and even pressure. Filling should occur from bottom upward to prevent voids. To facilitate placement of OPTICORE Z place a matrix band around the prepared tooth. (see diagram 8)

OPTICORE Z may be contoured by using a composite instrument. A flat-ended interproximal carver is recommended.

Place OPTICORE Z directly into the preparation and allow the system to self cure for 3:30 minutes. After that the material should be light-cured (40 seconds). With this technique an optimum of physical properties will be obtained.

An explorer can be used to test that the OPTICORE Z is completely set. Remove the matrix not earlier than the material has set.

Final crown preparation on the OPTICORE Z can be carried out by using crown preparation burs. (see diagram 9)

### Additional notes

- Do not use any resin to adjust viscosity of composite restorative material.
- Do not use in combination with substances containing eugenol because eugenol can impair the hardening of the composite and causes discoloration. Do neither store the composite material in proximity of eugenol containing products, nor let the composite allow to come into contact with materials containing eugenol.
- Contact of resin pastes with skin and gingival tissue should be avoided, especially by anyone having known resin allergies.

### Storage

Do not store above 20 °C (68 °F). Store unopened material in the refrigerator.

Opened cartridges have been used up within 3 months.

Do not use after expiry date.

### ATTENTION

This product is intended for professional dental use.

This product meets the requirements of **DIN EN ISO 4049**.

Selling of this product outside EU could be subject to limitations.

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