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GB OPTICORE PHOTO Light Cure Hybride Composite with Color Indicator

Rev2 – 3 May 2010

OPTICORE PHOTO is a visible light cure hybride composite for core-build-ups. OPTICORE PHOTO can be cured with 40 seconds by a light curing unit in a depth of cure of more than 4 mm. OPTICORE PHOTO is easy moldable and radiopaque. The filler content of OPTICORE PHOTO is about 83 %.

Preparation

1.1. Isolation

Use of a rubber dam to isolate the tooth is strongly recommended.

1.2. Cavity Preparation

Remove all existing old restorations and decay from the tooth. If necessary place any pins or posts (Opticore Post). Refer to directions of the selected post manufacturer. Prepare and clean the root canal with hydrogen peroxide (3%) or sodium hypochlorite solution (1-3%).

1.3. Pulp Protection

Cavity floor of deep excavations should be covered with a thin layer of hard setting calcium-hydroxide material.

1.4. Etching

Recommended is the total etch technique. Apply OPTIMA OPTIETCH onto the enamel

and dentine surfaces beginning with the enamel bevels. Leave the etching gel in place for **20 seconds**. On primary teeth and teeth high in fluoride, a **60 second** etch is recommended. Rinse for 20 seconds with water. Dry it in a water and oil-free airstream, but do not desiccate. A slightly wet dentine surface is important for the function of OPTICORE DUAL BOND ADHESIVE. The etched enamel bevel should have a chalky white appearance.

Etching Precaution:

It is essential, that etched areas are not contaminated by anything. If contamination occurs, re-etch, rinse with water and dry as above. Avoid etching gel contact with oral soft tissues, eyes and skin. If accidental contact occurs, flush immediately with copious amounts of water.

1.5. Application of OPTICORE PHOTO ADHESIVE

Apply OPTICORE DUAL BOND ADHESIVE generously with a brush onto the etched enamel and dentine surfaces and the fixed posts for 30 sec with agitation. Dry cautiously with oilfree air for 15 sec. Cure the coating with a dental halogen light unit for **20 seconds**.

Apply a second layer of OPTICORE DUAL BOND ADHESIVE generously with a brush

onto the etched enamel and dentine surfaces and the fixed posts for 30 sec with agitation. Dry cautiously with oilfree air for 15 sec. Cure the coating with a dental halogen light unit for **20 seconds**.

1.6. Application of OPTICORE PHOTO

OPTICORE PHOTO can now be applied for the build-up of the core. Using a non-metallic instrument is recommended.

1.7. Curing of OPTICORE PHOTO

Cure for **40 seconds** by using a light curing unit. OPTICORE PHOTO is filled to a high degree. Nevertheless to minimize the effects of polymerization shrinkage, curing in layers is recommended, if the the core-build-up is more than 4 mm deep. Leave for three minutes, proceed with trimming and finishing.

1.8. Storage

Do not store above 25 °C (78 °F). Cool storage prolongs shelf life.

Additional Notes

- Do not use any resin to adjust viscosity of composite restorative material.
- Contact of resin pastes with skin should be avoided, especially by anyone having known resin allergies.