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GB OPTICORE POST

Endodontic Tri-Glass fiber posts for chair side post and core build-ups

INSTRUCTION:

1 PREPARATION

NOTE: Use of rubber dam is recommended. Avoid using eugenol based temporaries.

NOTE: Cylindrical post with a size of 0.8 mm should only be used in cases of limited stress (anterior teeth or alongside other posts)

1. Post-canal check

After endodontic treatment, remove any temporary root filling and select an appropriate OPTICORE POST according canal size. Holding the post with tweezers seat it into the canal.

2. Root canal preparation

The post should fit into the canal without stress for a length of 6-7mm at least. The inserted post should stay firm without moving inside the root canal. Otherwise drill the root canal using the same drill size as the selected post.

NOTE: it is possible to use OPTICORE CANAL CLEANER as a drill lubricant.

3. Post shortening

After checking the post-canal fitting, determine the post length outside the canal according the estimated build-up size. Remove the post with tweezers and cut the post with diamond coated separating disc.

NOTE: Never use burs or cutters. It is suggested to clean the cut post to remove powder residues.

4. Canal cleaning and disinfection

Clean the canal using sodium hypochlorite and OPTICORE CANAL CLEANER. Leave this act for at least 25 seconds. Rinse thoroughly with sodium hypochlorite.

Dry gently with paper points leaving the surface slightly wet.

2 CEMENTATION

NOTE: OPTICORE POST is treated to present a micro rough surface and then covered by silane to increase bonding with cement. it is not necessary to apply any bond or primer onto OPTICORE POST

5. Adhesion

SELF ETCH - Apply OPTICORE ETCH BOND with a micro brush onto the slightly wet root canal walls for 15 seconds with agitation.

Repeat procedure 1 - 2 times. All the dentine surfaces must kept wet with the primer over the specified time.

Remove excess OPTICORE ETCH BOND with dry paper points. Dry the root canal cautiously with oil free air for 15 seconds to remove all volatile components and to disperse the adhesive to an even layer.

TOTAL ETCH - Etch enamel and dentine with OPTIETCH leaving in place for 15 seconds. Rinse thoroughly and dry gently with paper points leaving the surface slightly wet. Apply several layers of OPTICORE DUAL BOND mixed 1:1 with its dual activator. Dry the root canal gently with oil free air for 15 seconds to remove all volatile components and to disperse the adhesive all over the canal's walls.

6. Cementation

Insert the ENDO TIP into the canal until the apex seal is reached. Gently extrude the OPTICORE cement with a slow outward movement.

7. Post seating

Insert the post into the canal till it is completely seated. Light cure for 20 seconds in order to firm the post. Cementation will be complete after 3 minutes.

3 CORE BUILD-UP

8. Core build-up

Remove the ENDO TIP from the oral tip. Starting at crown base extrude OPTICORE always keeping the tip dip into the build-up material to avoid gaps entrapped into the restoration.

9. Finishing

OPTICORE will self cure in 3 minutes, anyway it is strongly recommended to light cure for 40 seconds in order to get optimal physical properties. Build-up finishing may be achieved with crown prep burs.

Additional notes

- Opticore posts are decontaminated by sodium hypochlorite during manufacturing.

- If necessary they can be cleaned and disinfected with isopropyl alcohol.

ATTENTION

This product is intended for professional dental use. Selling of this product outside UE could be subject to limitations.

International Dental Supply Srl

Via Valletta San Cristoforo 28/10 - 17100 Savona (SV) - ITALY

Tel +39 019-862080 - www.ids dental.it - email: info@ids dental.it