## Istruzioni d'uso 10 ml e (12 g):

### Descrizione

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## Instruction for use 10 ml e (12 g):

### Description

OPTIETCH EXTREME etching gel is a 37% phosphoric acid gel with optimum thixotropic flow properties. It is used for etching dental substances during enamel etching procedures and total-etch techniques. If the safety notes and instructions on use are observed, the product has no side-effects according to the current state of knowledge. Underfilling materials can be applied by the etching gel with a longer application time.

- Only for dental use!
- Do not use after the use-by date has expired.
- Close the syringe immediately after use.

### Preparation:

The area of the dental substance to be etched must be clearly marked.

- 1) The use of a rubber dam is recommended! Any contact with soft tissue (skin, oral mucosa, eyes etc.) is to be avoided. If etching gel gets into the soft tissue, rinse the area immediately with plenty of water and, if necessary, contact a doctor.
- 2) Before preparation, clean the tooth using a rubber cup with fluoride and fat-free polishing paste.
- 3) Remove cleaning residue carefully with water and clean and dry the surface of the tooth with oil-free compressed air.
- 4) Prepare the cavity accordingly, bevel the melted edges to achieve permanent enamel junctions. Rinse the cavity thoroughly with water and dry with oil-free compressed air. Avoid overdrying!
- 5) Protect dental pulp in deep cavities. (e.g. cover with a thin layer of suitable underfilling material). The cavity should be dry and free of contamination, whereby the overdrying of the dentine is to be prevented.

### Etching process

- 6) Apply the etching gel to the melted edges of the cavity and, where necessary, to the dentine.  
Recommended application time:  
15 - 30 sec. on dental enamel  
10 - 15 sec. on dentine  
50 - 80 sec. on milk teeth and  
Teeth with high fluoride content
- 7) Rinse for at least 20 sec. with water.
- 8) Dry with oil and water-free compressed air. The etched surface is not allowed to be touched until the treatment is started, or the surface has been contaminated by saliva. However, should contamination still occur after the etching and drying, the area must be etched again briefly (approx. 10 sec.).  
The etched enamel should look chalky white!  
After the etching, you can continue with the application of adhesives etc. in accordance with the manufacturer's instructions for use.

### Safety Instruction:

H314 Causes severe skin burns and eye damage.

Wear protective gloves/eye protection.

IN THE EVENT OF CONTACT WITH THE EYES: Rinse carefully for a long time with water.

Remove any contact lenses if possible. Continue to rinse.

Call a POISON CENTER or a doctor/physician immediately.

IF SWALLOWED: rinse the mouth. DO NOT induce vomiting.

IN THE EVENT OF CONTACT WITH THE SKIN: Wash with a lot of soap and water. Contains orthophosphoric acid.

ANGER

5°C  25°C

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IDS SpA - Via San Cristoforo 28/10 - 17100 Savona (SV) - ITALY  
Tel +39 019 862080 - [www.idsdental.it](http://www.idsdental.it) - email: [info@idsdental.it](mailto:info@idsdental.it)

Manufacturer: PRISMAN GmbH  
Otto-Hahn-Ring 6-18 - 64653 Lorsch  
+49(0)6251866980-00

CE 0124