

Blister da 2 frese

Blister da 3 frese

Blister da 5 frese

Depth Cutter 0.5mm	
Shape	834
Size (1/10mm)	021
Lenght (mm)	6.0
REF	59 02 01 M

Turbo					
Shape	368T	863T	850T	850T	850KRT
Size (1/10mm)	023	014	014	018	018
Lenght (mm)	5.0	8.0	10.0	10.0	9.0
REF	59 01 01 T	59 01 02 T	59 01 03 T	59 01 04 T	59 01 05 T

	Occlusal Reduction/Finishing				Wheels	Needle		
	Shape							
Shape	379	368	368	811	909	889	859	852
Size (1/10mm)	023	018	023	033	018	008	014	021
Lenght (mm)	4.2	4.5	5.0	4.0	0.5	4.0	10.0	10.0
XF Extra Fine			60 830EFFG023 \$17				60 859EFFG014 \$41	
F Fine 100µ	60 833FFG023 \$22	60 830FFG018 \$18	60 830FFG023 \$19			60 890FFG008 \$61	60 859FFG014 \$42	60 859FFG021 \$43
M Medium 130µ				60 811FG033 \$14		60 890FG008 \$62	60 859FG014 \$44	
C Coarse 145µ	60 833GFG023 \$23	60 830GFG018 \$20	60 830GFG023 \$11	60 811GFG033 \$15	60 815FG018 \$16		60 859GFG014 \$45	60 859GFG021 \$46
SC Super Coarse 180µ								

	Flame			Curettage	Round		
	Shape						
Shape	862	863	863	886	801	801	801
Size (1/10mm)	012	012	016	016	010	014	023
Lenght (mm)	8.0	10.0	10.0	10.0	1.0	1.4	2.3
XF Extra Fine	60 862EFFG012 \$47	60 863EFFG012 \$51				60 801EFFG014 \$1	
F Fine 100µ	60 862FFG012 \$48	60 863FFG012 \$52	60 863FFG016 \$53				60 801FFG023 \$2
M Medium 130µ	60 862FG012 \$49	60 863FG012 \$54		60 886FG016 \$59	60 801FG010 \$3	60 801FG014 \$4	
C Coarse 145µ	60 862GFG012 \$50	60 863GFG012 \$55	60 863GFG016 \$56	60 886GFG016 \$60	60 801GFG010 \$5	60 801GFG014 \$6	60 801GFG023 \$7
SC Super Coarse 180µ							60 801SGFG023 \$8

	Inverted Cone	Double Cone	Flat End Cylinder		Pear	Round End Cylinder
	Shape					
Shape	805	806	835	835	830	880
Size (1/10mm)	018	012	010	008	014	012
Lenght (mm)	2.3	3.0	3.5	3.0	5.0	8.0
XF Extra Fine						
F Fine 100µ						
M Medium 130µ		60 806FG012 \$10			60 808RLF014 \$11	
C Coarse 145µ	60 805GFG018 \$9		60 835GFG010 \$25	60 835GFG008 \$24	60 808RLGFG014 \$12	60 881GFG012 \$57
SC Super Coarse 180µ					60 808RLSGFG014 \$13	60 881SGFG012 \$58

	Round End Taper				Flat End Taper	
	Shape					
Shape	850	855	850	856	848	848
Size (1/10mm)	012	012	018	025	016	018
Lenght (mm)	10	6	10	7	10	10
XF Extra Fine	60 852EFFG012 \$33					
F Fine 100µ	60 852FFG012 \$34		60 852FFG018 \$35			
M Medium 130µ	60 852FG012 \$36		60 852FG018 \$37			60 848FG018 \$26
C Coarse 145µ	60 852GFG012 \$38	60 849GFG012 \$30	60 852GFG018 \$39	60 849GFG025 \$31	60 848GFG016 \$27	60 848GFG018 \$28
SC Super Coarse 180µ			60 852SGFG018 \$40	60 849SGFG025 \$32		60 848SGFG018 \$29

Blister da 5 frese

	Crown cutter					
	Shape					
Shape	CB21	CB31	CB31R	CB5CK	CB4MCL	CB35C
Size (1/10mm)	010	010	012	012	012	012
Lenght (mm)	4.1	4.1	4.1	4.0	4.0	8.0
FG (Gambo corto)	60 CB21R314010 \$63	60 CB31314010 \$64	60 CB31R314012 \$65	60 CB5CK314012 \$68	60 CB4MCL314012 \$67	60 CB35C314012 \$66
FG XL (Gambo lungo)				60 CB5CK316012 \$69		