

Dual phase teeth whitening system



COD. OC-22-0447-01
Code Final Product: 430447

ITEM	U.M.	LIMIT	TEST METHOD
VISUAL APPEARANCE	-	Pink pastel	Comparison with standards
PROFUM/FLAVOR	-	PINKOPOLIS - EAL52297/00 (menthol and red fruits)	Comparison with standards
UNIFORMITY	-	Homogeneous	Spread out on a slide
VISCOSITY TO	mPa.s	350.000 +/- 100.000	Internal method ANF001 gir7 1,5rpm 25°C
pH sol. 10%	-	7,50 +/- 0,50	UNI24003
DENISTY	g/ml	1,20 - 1,30	Internal method ANF013- ANF017
IONE FLUORIDE	ppm	920 - 1020	Method of external laboratory
TOTAL BACTERIAL COUNT	ufc/g	<10	UNI EN ISO 21149:2009
YEASTS AND MOULDS COUNT	ufc/g	<10	UNI EN ISO 16212:2011

PACKAGING	MATERIALS	ACCESSORIES
Tube 75 ml with four color printing + withe cap	Tube: ABL Cap: PP	Warranty seal

Ingredients:

Aqua, Hydrogenated starch hydrolysate, Hydrated silica, Glycerin, Potassium nitrate, Xylitol, Aroma, Sodium lauroyl sarcosinate, Cellulose gum, Sodium monofluorophosphate, Salvia officinalis leaf extract, Juniperus communis fruit extract, Calcium glycerophosphate, Hydroxyapatite, Menthol, Menthyl lactate, Sodium bicarbonate, Kaolin, Linseed acid, Sodium saccharin, PPG-26-buteth-26, PEG-40 hydrogenated castor oil, Sodium methylparaben, Potassium sorbate, Sodium benzoate, CI 77891, CI 16255.
Contains Sodium Fluoride: 1000ppm F-

RDA (Relative Dentin Abrasion):

CARE RDA 53

QC MANAGER SIGNATURE:

ISSUE: Quality check	VERIFY: Quality Ass. Manager	APPROVAL: Technical director

Dual phase teeth whitening system



COD. OC-22-0448-01
Code Final Product: 430446

ITEM	U.M.	LIMIT	TEST METHOD
VISUAL APPEARANCE	-	Light blu	Comparison with standards
PROFUM/FLAVOR	-	OUT OF THE BLUE EAL52301/00 (mentholated-fresh)	Comparison with standards
UNIFORMITY	-	Homogeneous	Spread out on a slide
VISCOSITY TO	mPa.s	400,0000 +/- 100000	Internal method ANF001 gir7 1,5rpm 25°C
pH sol. 10%	-	7,60 +/- 0,50	UNI24003
DENISTY	g/ml	1,20 - 1,30	Internal method ANF013- ANF017
IONE FLUORIDE	ppm	920 - 1020	Method of external laboratory
TOTAL BACTERIAL COUNT	ufc/g	< 10	UNI EN ISO 21149:2009
YEASTS AND MOULDS COUNT	ufc/g	< 10	UNI EN ISO 16212:2011

PACKAGING	MATERIALS	ACCESSORIES
Tube 75 ml with four color printing + with the cap	Tube: ABL Cap: PP	Warranty seal

Ingredients:

Aqua, Hydrogenated starch hydrolysate, Hydrated silica, Glycerin, Potassium nitrate, Tetrasodium pyrophosphate, Sodium lauryl sulfate, Aroma, Cellulose gum, Sodium monofluorophosphate, Zinc citrate, Menthol, Sodium saccharin, Papain, Sodium methylparaben, Disodium EDTA, Hydrogenated palm oil, PEG-200, Urea peroxide, Lecithin, CI 42090.

Contains Sodium Fluoride: 1000ppm F-

RDA (Relative Dentin Abrasion):

ACTIVE RDA 70

QC MANAGER SIGNATURE:

ISSUE: Quality check	VERIFY: Quality Ass. Manager	APPROVAL: Technical director